


## New Skin Tests at Intertek

In Compliance with REACH Directive, all cosmetics and chemicals will be tested by *in vitro* methods.

INTERTEK gives you the best way to analyze your products with Episkin® and Skinethic® human reconstructed epidermidis.

### Various Skin Models



- Corneal
- Oesophageal / Gingival / Vaginal
- Colored and non-colored skin...

### Applications

- Skin Irritation Assay (ECVAM validated protocol )
- Ocular Irritation Assay (pre-validated by ECVAM)
- Skin Corrosion testing (under TG 431 directive)
- Phototoxicity Testing of Raw Materials and Finished Products
- Bacterial Adhesion
- Biocides / Desinfectants Screening



### Contact Us

Campus industriel  
« Le Grand Chalon en Bourgogne »  
Route de Demigny  
71102 CHALON SUR SAONE Cedex  
Phone : 33-(0)3-85-99-12-76  
Fax : 33-(0)3-85-99-12-88  
[www.intertek-cb.com](http://www.intertek-cb.com)  
email : [chalon@intertek.com](mailto:chalon@intertek.com)

